



BY SAKURION

# CRYOCORE

NON-INVASIVE BODY SCULPTING

"CRYOLIPOLYSIS & BODY CONTOURING  
THE MODERN SOLUTION FOR SCULPTING YOUR IDEAL BODY"

# INTRODUCTION CRYOCORE MACHINE



NON-INVASIVE FAT FREEZING AND BODY SCULPTING

**Welcome** to the future of non-invasive fat reduction—our cutting-edge **Cryo-core** machine is designed to redefine body contouring with precision, safety, and effectiveness. Using advanced cooling technology, this innovative device targets and eliminates stubborn fat cells without surgery, downtime, or discomfort.

**Cryo-core** machine is the ideal solution for clients seeking noticeable results and a sculpted physique. Engineered with user-friendly features and tailored settings, it ensures a seamless experience for both practitioners and patients.

---

Cryolipolysis combined with electromagnetic technology freezes fat, reduces thickness, stimulates muscle growth, burns excess fat, and effectively shapes the body.



# sakurion®

[sakurion.com](http://sakurion.com)



# HOW DOES CRYOCORE MACHINE WORK ?

## Vacuum Massage\*

Prolonged suction during the treatment process causes discomfort and increases the risk of frostbite. Our patented vacuum massage\* technology solves this challenge in the beauty industry, ensuring safe treatments while providing customers with a more comfortable and relaxed experience.



1<sup>ST</sup> STEP:  
Pre-Heating



2<sup>ND</sup> STEP:  
360° Cooling



3<sup>RD</sup> STEP:  
Vacuum+Massage Mode



4<sup>TH</sup> STEP:  
Cryo Treatment


## TREATMENT APPLICATION




You may struggle with stubborn fat even after jogging and dieting.



The Cryo-Core suctions your skin and subcutaneous layer.




Cooling energy then freezes fat deposit to induce apoptosis.



Crystallized fat cells die off through the body's natural metabolism.



Remaining fat in the target area becomes leaner and less dense.



Combine with a healthy regimen to achieve a slimmer and more confident look.

UNIQUE PATENTED TECHNOLOGY ENSURES BEST RESULTS, COMFORT, AND SAFETY

# WHAT IS CRYOCORE ?


**CryoCore** combines advanced 360° Cryolipolysis and high-intensity focused electromagnetic technology for complete body contouring. The Cryolipolysis system naturally freezes and eliminates fat cells, targeting specific areas without damaging surrounding tissue. At the same time, high-intensity focused electromagnetic technology induces powerful muscle contractions, strengthening and toning muscles. With its 360° cooling technology and multi-applicators, CryoCore adapts to every curve, offering precise, non-invasive fat reduction and muscle sculpting.

## PERSONALIZED HANDPEICE

Not only 4 Handpeices can work at same time, User can also adjust cooling and suction levels independently to provide customized sessions for their Patients.

## FASTER TO TARGET POINT

The cooling system can reach its target temperature faster than previous and, maximizing time efficiency and convenience for both patients and practitioners.



**Applicators  
working  
in unison**

Multiple combinations  
in just one session.

## EFFECTIVE FAT REDUCTION

Effective fat freezing and muscle building with a combination of advanced cryolipolysis and high-intensity focused electromagnetic technology, enhancing muscle tone, reducing fat, and improving body contour for optimal results.

# FAST AND EASY APPLICATOR EXCHANGE

With a simple press-and-release mechanism, the **CryoCore** machine offers maximum convenience, making cooling cup changes effortless. Detaching cables or turning off the system is unnecessary, even during treatments.



## 8 INNOVATIVE DESIGNS TO REDUCE CELLULITE SAFELY AND MORE COMFORTABLY

### 7 IN 1 MULTI-SIZE APPLICATOR

Ergonomic design, skin-friendly fit the body arc Meet the whole body fat reduction needs from the face to the legs

### MEDICAL FLEXIBLE LSR SILICONE CUPS

Whole shaped molding 360° cups, delicate and soft, comprehensive care for the skin



### • PATENTED MASSAGE MODE

Effectively relieve the tension of the superficial muscles and speed up the blood circulation.

### • INNOVATIVE CONSTANT TEMPERATURE TECHNOLOGY

Independent research and development of constant temperature technology ensure precise temperature control and intelligent monitoring of temperature changes throughout the treatment.

### • UPGRADE VACUUM SYSTEM, THE SUCTION IS INCREASED BY 20%

A powerful negative pressure air pump allows for the suction force to be freely adjusted, enabling customized treatments.

### • SELF-CLEANING


Equipped with an intelligent self-cleaning function, the system offers one-button automatic cleaning after treatment, saving time and effort.

### • 4 APPLICATORS COULD WORK AT THE SAME TIME OR INDIVIDUALLY

Each handle has an independent engine, capable of operating individually or in combinations of 2, 3, or 4 simultaneously, delivering maximum power.

### INNOVATIVE "PREHEAT" TECHNOLOGY

Intelligent preheating system,  
Heating First and then cooling,  
Soothing the body,  
Enhancing blood circulation,  
Safe and Comfortable.



### 7 IN 1 MULTI-SIZE APPLICATOR

Ergonomic design, skin-friendly fit the body arc meet the whole body fat reduction needs from the face to the legs

### Medical flexible LSR silicone cups

Whole shaped molding 360 cups dedicated and soft, comprehensive care for the skin.

### 360 Cooling

Better Cooling effect and less treatment time.

### Liquid Silicone

Soft and tender touch like a baby skin.

360°

### Dual Temperature

Exclusive patent assures safety and efficacy.



### POWERFUL EMS

High-Intensity Focused Electromagnetic (HIFEM) technology uses electromagnetic induction to stimulate motor neurons, causing deep muscle contractions. This process enhances muscle density, reduces subcutaneous fat, and improves body contour in treated areas.

### TWO TYPES OF ERGONOMIC APPLICATORS


The patented curved applicator fits perfectly in hard-to-reach areas such as limbs, the waistline, and buttock contours.

### TREAT 4 AREAS SIMULTANEOUSLY IN ONLY 40 MINUTES

The **CryoCore** is equipped with an expandable range of applicator, making it easier than ever to achieve significant results on the most challenging of areas.

### UPGRADED TECHNOLOGY FOR ENHANCED RESULTS

The upgraded 360° surround cooling of the CryoCore delivers cooling energy evenly throughout the entire applicator, reaching the designated treatment temperature faster and increasing efficiency by up to 25%.



360° Surround cooling technology



# Technical specifications

## BASIC OPERATION

Technologies	high-intensity focused electromagnetic technology and Cryolipolysis
Handpiece	4 Handles : 2 EMS + 2 CRYO
User Interface	Intelligent Operation System
Display	15.6" HD Screen
Working Mode	4 Handles can work at same time or work independently.

## CRYO HANDLES TECHNICAL CHARACTERISTICS:

Maximum freezing temperature: -11°C, except for the Max Plus and Art-T models, where the maximum temperature is -9°C.

Preheating mode

Massage mode

Automatical self-cleaning

Emergency stop handle

## EMS HANDLES TECHNICAL CHARACTERISTICS:

Handle type	Flat handle or Curved handle
Treatment areas	abdomen, buttocks, arms, thighs, calves
Treatment time	0-60 minutes
Frequency	1-150 Hz
Tesla	up to 2.9T for each handle

**Voltage** 220V 50Hz

Working Temperature 5°C-40°C

**Treatment Area** Double Chin, Axilla, Knees, Arms, Inner Thighs, Outer Thighs, Abdomen, Flanks, Love Handles, Bra Line, Below Buttocks.



**sakurion®**

[www.sakurion.com](http://www.sakurion.com)

+32 2 315 45 43

[Info@sakurion.com](mailto:Info@sakurion.com)

Brand Whitlocklaan 154 1200 Sint-Lambrechts-Woluwe, Belgium

